

Molecular Diagnostics: Ahead of the curve with UNIFlow™

“I found the UNIFlow system to be one of the best Molecular Dx-SPECIFIC systems I have ever seen.”

CAP Mol. Pathology Lab Accreditation Inspector



For molecular diagnostics labs, where change is the only constant, UNICConnect delivers a solution that can keep pace. The core bio-molecular software technology behind UNIFlow makes it the perfect solution for Molecular diagnostic labs. UNICConnect provides labs that are on the cutting edge of developing science with a system that can keep up and keep them in compliance. If your lab is currently tracking information either in spreadsheets or in outdated systems that are difficult to manage, UNIFlow from UNICConnect will help you eliminate errors, increase productivity and significantly improve the quality of your results.

- ☺ Track samples from the clinic to the final report, and every step in between
- ☺ Flexibility to keep up with continuous process improvement
- ☺ Notifications to assure adequate inventory and personnel to run the assays
- ☺ Qualify every resource: Reagents, Workers, Instruments, SOPs
- ☺ Produce both internal and external reports
- ☺ Convenient access to get straightforward answers from data stored in multiple, complex databases

UNIFlow Customers and Compliance

UNIFlow customers consistently pass quality/accreditation inspections and audits with impressive results. UNIFlow evolves with the lab to stay current and meet all future requirements. Current customers are compliant with the following:

- | | | | |
|-------|--------|-------------|------------------|
| ☺ GMP | ☺ FDA | ☺ ISO 17025 | ☺ 21 CFR Part 11 |
| ☺ CAP | ☺ CLIA | ☺ GLP | ☺ AABB |

UNIFlow is a Hybrid System Made From Pre-Built and Custom Modules Designed for MDx Labs

- Pre-built: Accession, Extraction, Quantify, PCR, Analysis, Results
- Custom: Internal and External Reports, Unique Steps, Input Forms
- Use familiar names, labels and titles
- Input screens are intentionally kept simple and straightforward

Bar code background prompts technician to scan rather than hand key information

Other step groups or modules are always just a click away

Queues drive the process and shows where samples are currently

Contact UNIConnect to see a live system demo via the Web!